

# Applications of Deep Learning for High-Throughput Imaging

Gianluca Pegoraro  
CCR High-Throughput Imaging Facility  
Jun 15 2021

# High-Throughput Imaging (HTI)

PI


Experimental  
perturbation

Imaging-based  
cellular assay

Phenotypic  
change

HiTIF


Automated  
liquid handling



High-throughput  
microscopy









High-content  
image analysis



Up to:

- $10^4$  Wells
- $10^4$  Cells/Well
- $10^2$  Feat./Cell

# High-Throughput Acquisition and Analysis




Results

$$2D \text{ images/day} = n * m * \lambda * z * t \approx \text{up to } 2 * 10^5$$


# Deep Learning for Nucleus Segmentation

- Accurate Detection:
  - 90%-95% accuracy
- Practical:
  - Trainable with ~ 10 FOVs (~500 - 1,000 objects)
  - Fast inference (~ 1s/FOV)
- Robust and Generic:
  - Different cell types
  - Different magnifications
  - Different confluency


# Semi-automated GT Label Generation




# Image Augmentation and Bootstrapping



# Feature Pyramid Networks (FPN)-Watershed (WS)




# Pipeline for Training and Testing DL Models






# DL Models Trained on MCF10A Images (1)




# DL Models Trained on MCF10A Images (2)




# F1 Score to Measure Inference Performance

$$F1(t) = TP(t) / (TP(t) + (FP(t) + FN(t)) / 2)$$



IoU(Threshold) = 0.50




Ground Truth





Prediction





# Inference Performance of Baseline DL Models




# Transfer Learning Improves MRCNN Performance




# Image Augmentation is not Required



# Final DL Models (1)




# Final DL Models (2)






# Final Models Performance




# Summary 1)

- Semi-automated computational pipeline for DL models training/testing
- Transfer learning can improve performance by using networks weights obtained from training on everyday objects
- Training vs. out of the box: it depends...
- Other DL applications: classification, denoising, inpainting

# Future Areas of Improvement for DL in Bioimaging



# CEM500K




# Imjoy

ImJoy



Embedded in Notebooks



Corrected Annotations

Predicted Labels

Interactive Model Trainer

Local/Remote Server




binder  
colab





TensorFlow


PYTORCH




# CellPose: 2D Segmentation



# CellPose Works "Out of the box"




# CellPose: 3D Segmentation





# Better Tools to Serve Models: Deep Cell Kiosk

## DeepCell Deployment Kiosk Architecture



## Summary 2)

- Rapid improvements in making DL more accessible for biologists, larger curated datasets, better model architectures, higher-throughput at inference
- Biologists should pair up with ML/DL experts

# Acknowledgements

## HiTIF

P. Gudla (f)

L. Ozbun

K. Lee (f)

A. Keikhosravi

## FNL

G. Zaki

J. Kim

## Collaborators

S. Shachar (f)

T. Misteli (CCR)

M. Gadkari

L. Franco (NIAMS)

J. Sung

I. Fraser (NIAD)

NIH HPC Team

# Contact Information



*Laboratory of Receptor Biology and Gene Expression*

*Center for Cancer Research/NCI*

Gianluca Pegoraro

[gianluca.pegoraro@nih.gov](mailto:gianluca.pegoraro@nih.gov)